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1. Burgeson RE. New collagens, new concepts. *Ann Rev Cell Biol* 1988;4:551-577
2. Fledelius C, Johnsen A, Cloos P, Bonde M, Qvist P. Identification of a β -isomerized aspartyl residue within the c-terminal telopeptide $\alpha 1$ chain of type I collagen. Possible relation to aging of bone. *J Bone Miner Res* 1996;11(Suppl.1) Abstract No. 113.
3. Bonde M, Qvist P, Fledelius C, Riis BJ, Christiansen C. Applications of an Enzyme Immunoassay for a New Marker of Bone Resorption (CrossLaps): Follow-up on Hormone Replacement Therapy and Osteoporosis Risk Assessment. *J Clin Endocrinol Metab* 1995;80:864-868.
4. Ravn P, Clemmesen B, Riis BJ, Christiansen C. The Effect on Bone Mass and Bone Markers of Different Doses of Ibandronate: A New Bisphosphonate for Prevention and Treatment of Postmenopausal Osteoporosis. A 1-year, Randomized, Double-Blind, Placebo-Controlled Dose-Finding Study. *Bone* 1996;19(5):527-533.
5. Garnero P, Sornay-Rendu E, Duboeuf F, Delmas F. Markers of Bone Turnover Predict Postmenopausal Forearm Bone Loss Over 4 Years: The OFELY Study. *J Bone Miner Res* 1999;14:1614-1621.
6. Ganero P, Borel O, Delmas PD. Evaluation of a Fully Automated Serum Assay for C-Terminal Cross-Linking Telopeptide of Type I Collagen in Osteoporosis. *Clin Chem* 2001;47(4):694-702.

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