

Guidelines for the use of blood warming devices (AABB)

Reference: 2002 AABB





Introduction

- Hypothermia
 - Induced by rapid, large-volume transfusion of refrigerated blood components
 - A potential source of serious complications
 - Especially in pediatric and elderly
- Prevent systemic hypothermia
 - Warming of blood component and other IV fluids during massive transfusion




Blood warming technologies

- Countercurrent heat exchange
- Dry heat
- Thermostatically controlled waterbath
- In-line microwave



- Countercurrent heat exchange
 - Have a multiple lumen tubing assembly
 - Water circulates the length of the tubing to the patient
 - Preventing conductive heat loss before infusion
 - Have a fluid delivery temperature of 35~41 C with rapid flow rates and low priming volumes




Reference : Canadian journal of anaesthesia, 1995






- Dry heat
 - Previous comparative studies indicated less efficiency
 - Recent improvements and design changes
 - More convenient to use
 - No water chambers are required
 - Some do not require special transfusion sets.



Reference : Canadian journal of anaesthesia, 1995





- Thermostatically controlled waterbath
 - Immersing the entire unit of blood into the heated water
 - This is older technology and has several safety concerns
 - To ensure even warming, blood must be agitated during the warming process
 - The unit is hung for transfusion it begin to cool rapidly
 - It is possible to inadvertently return the unit to general inventory should the unit not be needed.



- In-line microwave
 - Single-patient, disposable-tubing assemblies that fit into the warming chamber (where the blood passes through, is heated to approximately 37 C)
 - Affects the warming efficiency
 - The initial temperature of the blood component
 - The flow rate through the heating chamber



- The blood is heated while passing through a small disposable cartridge inserted into the microwave cavity.
- The microwave energy emitted in the cavity is strictly controlled
- The temperature of the fluid is measured and displayed throughout the heating process



- Blood warmers
 - Must have an **exact temperature control** to prevent hemolysis caused by excessive heat
 - An **alarm system** and a **visible temperature display** are required
 - In order to comply with AABB, CAP, FDA regulation, a blood warmer must be **properly maintained** and either be **FDA-approved**



- As with all equipment used in the blood bank and transfusion process
 - Must be procedures for equipment operation
 - Validation
 - Routine preventive maintenance
 - Quality control(training and competency of staff)
 - Properly documented installation, validation, and maintenance records
 - The blood bank medical director should be actively involved in establishing these procedures



Advantages of blood warming devices

- To prevent the complication caused by hypothermia
 - Improving survival rate and patient outcomes (including decreased length of hospitalization)
- Hypothermia
 - Defined
 - Reduction in core body temperature below 35C



Hypothermia

- Side effects
 - Cardiac arrhythmia
 - Hemostasis abnormalities from impaired platelet function and slowed enzymatic reactions in the coagulation cascade
 - Vasoconstriction
 - Dehydration
 - Lack of oxygen to tissue
 - Increased red cell release of potassium
 - Citrate toxicity(with blood component transfusion)



- Side effects
 - The metabolism of drugs is also impaired(prolong)
 - Immune function impaired (reduced oxygen in tissue impairs oxidative killing of bacteria by neutrophils)
 - May promote surgical-wound infection (reduces the deposition of collagen during wound healing)



- If the core body temperature reaches 25~30 C
 - Cardiac output decreases
 - Ventricular irritability increase
 - Possibly leading to ventricular fibrillation
 - Cardiac arrest



Disadvantages and complications

- Risk of hemolysis
- Risk of sepsis
- Decreased Infusion Rate
- Inefficient Heat Transfer
- Expense



- Risk of hemolysis

- The red cell membrane become less pliable at higher temperature and is more susceptible to the effects during rapid transfusion
- Thermal injury to red cells may occur if is not used according to the manufacturer's directions
- If the flow within the warmer is stopped temporarily, the temperature may rise above 37C, resulting in slight hemolysis
- Out of order



- Risk of sepsis
 - Due to bacterial growth at warmed temperature is increased (especially if the blood is kept warm for a prolonged period of time)
 - All disposables should be inspected before use and installed according to the manufacturer's directions (damaged disposables could leak contaminated water)



- Risk of sepsis(2)
 - Warm waterbaths
 - Have the potential risk of harboring pathogens that could contaminate the ports or injection sites of the blood bags
 - Water may be splashed into the sterile field during insertion of the transfusion disposables
 - Water-tank decontamination procedure is essential



- Decreased Infusion Rate
 - The long, narrow tubing increases the resistance to flow and reduces flow rate.
 - In trauma situation, this decrease in infusion time may be detrimental
 - The selection of equipment designed for rapid infusions will offset this drawback



- Inefficient Heat Transfer
 - Limited by the flow rate of the infusion
 - Limited by the heating technology of the device in use
 - Most do not deliver normothermic (37C) fluids at very rapid flow rate(only 33~36C)



- Expense
 - The cost of the equipment
 - The disposables
 - Staff training
 - Quality control
 - Preventive maintenance



Indication for use

- Massive transfusion
- Trauma situations in which core-rewarming measures are indicated
- Administration rate $>50\text{mL/min}$ for 30 minutes or more (Adult)
- Administration rate $>15\text{mL/kg/hour}$ (Children)
- Exchange transfusion of a newborn
- Patient re-warming phase during cardiopulmonary bypass surgical procedures



- Using blood warmer should be based on
 - The size and clinical status of the patient
 - Blood loss
 - Replacement needs
 - Rate of transfusion
 - Available warming systems



- Other clinical situations where blood warming may be considered
 - Potent, high-titered, cold autoantibodies, reactive at body temperature and capable of binding complement, thus causing hemolytic anemia
 - Raynaud's syndrome or other cold-induced vasoactive effects
 - Neonatal and pediatric transfusions
 - Therapeutic plasma apheresis or red cell exchange procedures



Contraindications

- Blood warmers are not usually needed or provide no clinical advantage in situations such as :
 - Elective transfusions at conventional(slow) infusion rate
 - Most cold agglutinins encountered in pretransfusion testing
 - Patient experiencing shivering and discomfort due to the cold(Method to warm the patient, not the blood)
 - Administration of Platelets, Cryoprecipitate, or Granulocyte suspensions (Warming may render these products less effective)



Warnings

- There many warnings to consider when using blood warming equipment:
 - Do not use fluid warming equipment that is not FDA-approved for use with blood components, or is not internally validated
 - Do not use faulty blood warming equipment
 - Do not use transfusion sets not approved for use with the blood warming device
 - Do not warm blood components by placing on or near a radiator, heater, stove, patient-warming blanket, or in a conventional microwave oven or plasma thawer



Warnings(2)

- Do not allow the unit to sit at room temperature for prolonged periods to warm up
- Do not place blood components under running hot tap water or in an unmonitored or improvised warm waterbath
- Do not place pediatric syringe aliquots or transfusion tubing under warming lamps. Validated shielding should be used to prevent inadvertent and uncontrolled (over) warming when tubing has to pass through a warmed field



Warnings(3)

- Do not return blood components that have been warmed to inventory
- Even warming is imperative, the addition of heated saline as a diluent/admixture to red cell is difficult to standardize or control
- However, it does have the advantages of rapid heating, not requiring expensive equipment, and the rapid infusion of diluted blood



Administrative aspects

- Choosing an appropriate blood warming device should be based on
 - Equipment cost
 - The cost of disposables
 - The ease of use(Training concerns)
 - The disposables used should be latex-free
 - Noted in the transfusion record




Compliance Issues(AABB)

- Blood warming devices must have
 - A visible thermometer
 - A warning system to detect malfunction
 - The warming method should not cause hemolysis
- The standards regarding critical equipment and process control include:



- Defined processes for installation, calibration, validation, and documentation of problems encountered
- Equipment management
- Development of procedures and schedules for calibration and maintenance
- Documentation of period calibration and maintenance, timeliness of repairs, and post-repair calibration



Compliance Issues(CAP)

- TRM41500
 - An FDA-cleared blood warming system must be equipped with special feature to alert the user to propertransfusion conditions
 - The blood-warming system must be
 - FDA-cleared
 - Properly maintained
 - Special features(a visible thermometer and a audible alarm)
 - Does not damage to the blood component being warmed



- TRM41000
 - Transfusion service must have a written policy and procedure for blood administration
 - Transfusion personnel have been trained and are competent in the procedure



Quality Assurance

- Ownership
 - The primary difficult encountered by transfusion services
- Several different brands or types of blood warmers may be in use within the hospital
- Distributed through several departments such as surgery, emergency.....
- The users may include nurses, anesthesiologists, technicians.



- The transfusion service medical director responsibility
 - Ensuring the proper function of blood warmers
 - Determining procedures
 - Frequency of testing
 - Review and retention of records



- If the biomed department routinely performs only electrical safety checks:
 - Alarm activation
 - Effluent temperature/thermometer calibration
 - Cleaning and disinfecting of the water tank(If applicable)



- Quality control procedures and routine maintenance must be performed, documented and reviewed
- Any devices found to be out of control or faulty must be sequestered and promptly followed by corrective action
- Staff training and continued competency assessments must be performed and documented

Thank you for your attention

