




Cell Apoptosis/ Cell Cycle/Proliferation CBA 微珠免疫分析

Grace Lin 林秀鴻
Senior Product Specialist,
BD Biosciences
Grace.lin@bd.com


© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.

BD Pharmingen Total Solution




2 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.






Apoptosis

4 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company. 

Apoptosis


細胞凋亡，又稱為程式性細胞死亡 (Programmed Cell Death) 是機體主動的、高度有序地清除無用細胞的過程。



5 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Apoptosis




1. Phosphatidyl Serine exposed
2. Mitochondrial potential decreases
3. Caspases activated
4. DNA fragmentation

6 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




Apoptosis: Annexin V mechanism



7 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Annexin V/PI detect apoptosis



BD Pharmingen™ FITC Annexin V Apoptosis Detection Kit I (Cat.No. 556547) on Accuri C6

Live cells:
Annexin V Negative
PI Negative

Early Apoptotic:
Annexin V Positive
PI Negative

Late apoptotic/dead:
Annexin V Positive
PI Positive

8 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Protocol for Annexin V/PI

- Wash cells twice with cold PBS and then resuspend cells in 1X Binding Buffer at a concentration of 1×10^6 cells/ml
- Transfer 100ul of the solution to a 5ml culture tube
- Add 5ul of FITC Annexin V and 5ul PI
- Gently vortex the cells and incubate for 15min at RT in the dark
- Add 400ul of 1X Binding Buffer to each tube. Analyze by flow cytometry within 1hr.

Apoptosis kit:

556547	Annexin V FITC Apoptosis Detection Kit I
559763	Annexin V PE Apoptosis Detection Kit I

9 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




Tips:

- Optimized for suspension cells
- Compatible across most species
- Compatible for cell surface staining as long as the calcium concentration is maintained
- Using Binding buffer during the whole assay
- Check trypsin/EDTA for adherent cells
- Calculate compensation using single stained apoptotic cells
- Tricky to determine FITC-%PE compensation value sometime
- Negative group and biological control group
- PE Annexin V and APC Annexin V for GFP transfected cell lines

10 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Apoptosis




- 1. Phosphatidyl Serine exposed
- 2. Mitochondrial potential decreases
- 3. Caspases activated
- 4. DNA fragmentation


11 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




Apoptosis: Mitochondrial membrane potential



Thymocytes, untreated



Thymocytes, + dexamethasone



BD Pharmingen™ BD™ MitoScreen (JC-1)
(Cat.No. 551302) on Accuri C6


12 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Protocol for JC-1



BD Pharmingen™ BD™ MitoScreen (JC-1) contains 4 amber vials of lyophilized JC-1 and 60ml 10X assay buffer



13 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




Tips:

- Warm 10X Assay buffer at 37°C to completely dissolve salt crystals
- Keep diluted 1X Assay buffer at 37°C prior to use
- JC1 should be used as fresh as possible, avoid frozen and thaw repeatedly
- After dilute JC1 stock with 1X assay buffer, vortex working solution thoroughly
- After vortexing, JC1 aggregates may be present, but should not interfere with flow cytometric analysis.

14 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Apoptosis




- 1. Phosphatidyl Serine exposed
- 2. Mitochondrial potential decreases
- 3. Caspases activated
- 4. DNA fragmentation


15 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Apoptosis: Active Caspase-3



Nature Reviews | Microbiology



BD Pharmingen™ PE Active Caspase-3 Apoptosis Kit (Cat.No. 550914) on Accuri C6

550914	BD Pharmingen™ PE Active Caspase-3 Apoptosis Kit
560901	BD Pharmingen™ FITC Rabbit Anti- Active Caspase-3

16 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Protocol for active caspase-3 detection

- Wash cells twice with cold 1X PBS, then resuspend cells in BD Cytofix/Cytoperm solution at a concentration of 1×10^6 cells/0.5ml
- Incubate cells for 20min on ice
- Pellet cells, aspirate, and discard BD Cytofix/Cytoperm solution, wash twice with BD Perm/Wash buffer (1X) at a volume of 0.5ml buffer/ 1×10^6 cells at room temperature
- Resuspend each 1×10^6 cells in 100ul BD Perm/Wash buffer (1X) and 20ul antibody and incubate for 30min at room temperature
- Wash each test in 1.0ml BD Perm/Wash buffer (1X), then resuspend the test in 0.5ml BD Perm/Wash buffer (1X) and analyze by flow cytometry


550914 BD Pharmingen™ PE Active Caspase-3 Apoptosis Kit

560901 BD Pharmingen™ FITC Rabbit Anti- Active Caspase-3

17 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Apoptosis




1. Phosphatidyl Serine exposed
2. Mitochondrial potential decreases
3. Caspases activated
4. DNA fragmentation

18 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Apoptosis: DNA fragmentation



BD APO-DIRECT Kit Component

Part A:

- PI/RNase Staining Buffer
- Reaction Buffer
- Rinsing Buffer
- Wash Buffer


Part B:

- FITC-dUTP
- Negative Control Cells
- Positive Control Cells
- TdT Enzyme

19 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Apoptosis: DNA fragmentation



20 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Apoptosis product list

Cat.No.	Description
556547	Annexin V FITC Apoptosis Detection Kit I
559763	Annexin V PE Apoptosis Detection Kit I
551302	BD™ MitoScreen (JC-1)
550914	BD Pharmingen™ PE Active Caspase-3 Apoptosis Kit
560901	BD Pharmingen™ FITC Rabbit Anti- Active Caspase-3
556405	BD Pharmingen™ APO-BRDU™ Kit
556381	BD Pharmingen™ APO-DIRECT™ Kit
550378	BD Pharmingen™ Z-DEVD-FMK, Caspase-3 Inhibitor

21

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




Viability

22

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Viability staining improves results




Dye	Unfixed cells	Fixed cells	Detector	Laser
DAPI	✓	✗	BV421	UV/Violet
Via-Probe Green	✓	✗	FITC	Blue
PI	✓	✗	PE	Blue/YG
7-AAD	✓	✗	PerCP-Cy ⁵ .5	Blue/YG
DRAQ7™	✓	✗	APC	Red
Via-Probe Red	✓	✗	APC	
FVS450	✓	✓	BV421	Violet
FVS510	✓	✓	BV510	Violet
FVS575V	✓	✓	BV605	Violet
FVS520	✓	✓	FITC	Blue
FVS570	✓	✓	PE	Blue/YG
FVS620	✓	✓	PE-CF594	Blue/YG
FVS660	✓	✓	APC	Red
FVS700	✓	✓	AF700	Red
FVS780	✓	✓	APC-H7	Red

23 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



PI & 7-AAD: Cell membrane integrity



Live

Dead


無法針對固定過的細胞分別死活→解決方案:FVS染劑

24 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Fixable viability stain (FVS)-可固定死活染劑


- Amine-reactive membrane impermeable dyes useful for live/dead discrimination
- Labeled cells can be fixed and permeabilized, making the dyes compatible with multiple downstream applications 可與不同下游應用相容



25 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Fixable viability stain and apoptosis



26 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Protocol for FVS staining

- Harvest cells and wash once with azide- and protein-free 1X DPBS
- Resuspend cells at $1-10 \times 10^6$ cells/ml in 1X DPBS
- Add 1ul of Fixable Viability Stain 450 stock solution for each 1ml of cell suspension and vortex immediately
- Incubate the mixture for 10-15min at RT protected from light
- Wash cells twice with 2ml of Stain buffer
- Decant the supernatant and gently mix the cell pellet, resuspend the cell in Stain buffer
- Stain, fix and permeabilize cells as desired for down stream applications

27 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Tips for FVS staining

- Incubate mixture at 2-8°C for 20-30min for mouse cell applications
- Incubate mixture at 37 °C for 5-7min for BD Phosflow applications
- The reactivity of free dye is quenched by washing with buffer containing protein (FBS or BSA) prior to staining with fluorescent antibodies

28 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Viability product list

Description	Cat. No.	Laser
BD Pharmingen™ 7-AAD	559925	Blue
BD Pharmingen™ Propidium Iodide Staining Solution	556463	Blue
BD Horizon Fixable Viability Stain 520	564407	Blue
BD Horizon Fixable Viability Stain 660	564405	Red
BD Horizon Fixable Viability Stain 700	564997	Red
BD Horizon Fixable Viability Stain 780	565388	Red

29 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.





Cell cycle and Proliferation

30 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Cell cycle: DNA content




31 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



DNA probes

- PI-Propidium Iodide
- EB-Ethidium Bromide
- 7AAD-7-amino-actinomycin D
- AO-Acridine Orange
- Hoechst 33342
- DAPI- 4',6-diamidino-2-phenylindole
- DRAQ7
- Via-Probe Red



32 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Protocol for EtOH fixed PI/RNase staining

- Harvest cells and prepare single cell suspension in PBS
- Wash cell twice with PBS and resuspend at 1×10^6 cells/ml
- Aliquot 1ml cells in a 15ml PP V-bottomed tube and add 3ml cold absolute ethanol dropwise while vortexing
- Fix cells for at least 1 hour at 4°C
- Wash cells twice in PBS
- Add 0.5ml/ 10^6 cells of PI/RNase staining buffer and mix well. Incubate 15min at RT
- Store samples at 4°C until analyzed by flow cytometry

33 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




Protocol for Cycletest Plus DNA Kit

- Centrifuge the cell suspensions at 400g for 5min at RT
- Carefully decant all the supernatant
- Add 250ul of Solution A (trypsin buffer) to each tube and gently mix by tapping the tube by hand. DO NOT VORTEX
- Incubate for 10min at RT
- Add 200ul of Solution B (trypsin inhibitor and RNase buffer) to each tube and gently mix by tapping the tube by hand
- Incubate for 10min at RT
- Add 200ul of cold Solution C (PI stain solution) to each tube and gently mix tapping by hand
- Incubate for 10min in the dark on ice or in the refrigerator
- Filter the sample through 50um nylon mesh and then analyzed by flow cytometry


34 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Gating strategy for Cell cycle assay



1. FSC/SSC:
cell population
2. PI-A/PI-H(W):
double-let
discrimination
3. PI-A(H) Histogram:
cell cycle
distribution



35

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Tips:

- Fixation process using ethanol is critical for cell cycle assay
- Make sure discriminate doublelets
- Using **LOW** flow rate to decrease CV
- Record more than 20,000 cells in singlelets gate
- Using ModFit or other software to analyze cell cycle data


36

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



DNA Synthesis: BrdU


- BrdU is incorporated into the DNA of newly synthesized cells (S phase)
- Incorporated BrdU is stained with specific anti-BrdU antibodies
- Staining with a dye that binds total DNA is often coupled with BrdU



37 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Protocol for BrdU Kit




38 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Proliferation dye: CFSE

- CFSE (carboxyfluorescein succinimidyl ester) is non-fluorescent molecules that diffuse into cells. Once inside cells, they are hydrolyzed by intracellular non-specific esterases and covalently bind to cellular components to become fluorescent products



39 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Tips:

- Avoid staining in azide-, serum-, or protein-containing buffers during CFSE labeling
- Higher initial staining intensities can be achieved by increasing the initial CFSE concentration, but may lead to increase cell toxicity or cell death

40 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




Cell cycle and proliferation product list

Cat.No.	Description
340242	BD Cycletest™ Plus DNA Reagent Kit
550825	BD Pharmingen™ PI/RNase Staining Buffer
561980	BD Pharmingen™ Hoechst 33342 Solution
564907	BD Pharmingen™ DAPI Solution
565082	BD Horizon™ CFSE
557892	BD Pharmingen™ APC BrdU Flow Kit

41 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



DNA damage causes caspase independent Cell Cycle arrest




H2AX(Ser139) is a key to detect DNA double-stranded breaks (DSBs)

42 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Time course of γ H2AX fluorescence after UV irradiation



PNAS, June 2006, Vol 103, 9891-9896

43 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



γ H2AX products

Cat.No	Description
564720	BD Horizon™ BV421 Mouse Anti-H2AX (pS139)
562377	BD Pharmingen™ PE Mouse Anti-H2AX (pS139)
560447	BD Pharmingen™ Alexa Fluor® 647 Mouse anti-H2AX (pS139)
560445	BD Pharmingen™ Alexa Fluor® 488 Mouse anti-H2AX (pS139)
562253	BD Pharmingen™ Apoptosis, DNA Damage and Cell Proliferation Kit

44 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




Apoptosis, DNA Damage and Cell Proliferation Kit (Cat.No 562253) Components

Part	Components	Size
Part A	PerCP-Cy™5.5 Mouse Anti-BrdU	50T
	Alexa Fluor® 647 Mouse Anti-H2AX (pS139)	50T
	PE Mouse Anti-Cleaved PARP (Asp214) Antibody	50T
	BD Cytofix/Cytoperm™ Fixation/Permeabilization Solution	25ml
	BD Perm/Wash™ Buffer (10X)	25ml
	BD Cytofix/Cytoperm™ Plus Permeabilization Buffer	10ml
	DAPI	100ul
Part B	BudU(10mg/ml)	5mg
	DNase	300ul

45 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Apoptosis, DNA Damage and Cell Proliferation Kit (Cat.No 562253) protocols and data



46 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




Cell function assay example

47 © 2016 BD, BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Camptothecin induced Jurkat cell apoptosis



1X10⁶ Jurkat cells were treated with 8 μM camptothecin for 6 hours, harvest cells and wash twice with cold PBS, then perform assays according protocols provided by Cat.No. 556547, Cat.No. 551302, Cat.No. 550914 and Cat.No. 340242 (BD Biosciences). Data from BDB SH COE, acquired on BD FACSCanto II, analyzed by BD FACSDiva software and ModFit software.

48 © 2016 BD, BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Example: CPTH6 induced apoptosis in U937 cell on Accuri C6

Cancer Therapy: Preclinical

Clinical
Cancer
Research

CPTH6, a Thiazole Derivative, Induces Histone Hypoacetylation and Apoptosis in Human Leukemia Cells

Daniela Trisciuglio¹, Ylenia Ragazzoni¹, Andrea Pelesi², Marianna Desideri¹, Simone Carradori¹, Chiara Gabellini¹, Giovanna Maresca¹, Riccardo Nescatelli¹, Daniela Secci¹, Adriana Bolasco¹, Bruna Bizzanti¹, Chiara Cavallero¹, Igna D'Agnano¹, Patrizia Filetti¹, Lucia Ricci-Vitiani², Maria Giulia Rizzo², and Donatella Del Bufalo¹


Clin Cancer Res. 2012 Jan 15;18(2):475-86



49 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




CPTH6 is a specific Gcn5/pCAF inhibitor



50 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




CPTH6 reduces in vitro cell viability of U937 cells



51 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




CPTH6 decrease cell viability of NSCLC cells in vitro



52 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




CPTH6 induces cell-cycle perturbation in U937 cells



Flow cytometric analysis of DNA content:
Determination of fluorescent intensity of propidium iodide intercalated in DNA.


Example of profile:




53 © 2016 BD, BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



CPTH6 induces apoptosis in U937 cells




Flow cytometric TUNEL assay:
Apoptosis detection by a two-color assay for labeling DNA breaks and total cellular DNA content.




54 © 2016 BD, BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




Determination of apoptotic mechanisms induced by CPTH6 treatment in U937 cells




Western Blot analysis:
Detection of apoptotic markers



Flow cytometric assay detecting Reactive Oxygen Species production (DHE)




Extraction and Western Blot analysis:
Detection of an apoptotic marker in cytosol



55 © 2016 BD, BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




CPTH6 induces apoptosis in patient-derived stem-like LCSC136 cells



56 © 2016 BD, BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




CPTH6 induces specific alterations to the plasma membrane, hallmark of apoptotic cells



Flow cytometric AnnexinV binding assay:
Fluorescent Annexin V binding to altered plasma membrane of apoptotic cells.


Example of profile:



57 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.





CPTH6 induces Caspase-3 activation and cleavage of PARP in LCSC136



Flow cytometric active caspase 3 detection assay
Detection of the caspase 3 active form using an anti-active caspase-3 antibody.

Example of profile




Western Blot analysis:
Detection of cleaved PARP protein, an apoptotic marker

58 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




CPTH6 does NOT induce apoptosis in established human H1299 NSCLC cells



Flow cytometric AnnexinV binding assay:

Fluorescent Annexin V binding to altered plasma membrane of apoptotic cells.

Example of profile:




Western Blot analysis:

Detection of cleaved PARP, an apoptotic marker.

59 © 2016 BD, BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




CPTH6 induces DNA damage in established human H1299 NSCLC cells



Flow cytometric γ H2AX assay

Detection of H2AX phosphorylation, a marker for DNA damage, using a H2AX antibody.



Immunofluorescence analysis of H2AX phosphorylation

60 © 2016 BD, BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Summary

- Keep cells at good condition
- It is critical to determine IC50 value
- Since cell lines and drugs vary, concentration of drugs and time of treatment should be titrated
- Use multiple approaches to detect apoptosis

61 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



CBA Cytometric Bead Array

62 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Tools and Techniques for Cytokine Quantitation


Tool/Technology	BD Cytometric Bead Array (CBA)	ELISA	ELISPOT
Molecules detected	Secreted or intracellular	Secreted	Secreted (in situ)
Multiparameter	Yes	No	No
Single-cell/cell-subset information	No	No	Frequencies, no subset information
Antigen-specific	Yes	Yes	Yes
Post-assay viability	Yes	Yes	No
Protein quantitation	Yes	Yes	No
Instrumentation	Flow cytometer	Spectrophotometer	ELISPOT reader

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




ELISA

Direct



Sandwich





© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



ELISPOT


- 1 Capture Antibody**
Coat ImmunoSpot™ plate micro well with anti-cytokine capture antibody
- 2 Blocking**
Block unoccupied well sites with protein
- 3 Add Cells**
Incubate cells in well
- 4 Wash**
Cells are washed off

- 5 Detection Antibody**
Add Biotinylated anti-cytokine detection antibody
- 6 Enzyme-Avidin**
Add Avidin-HRP
- 7 Develop With Substrate**
Add substrate and monitor formation of colored spots

© 2016 B

property of Becton, Dickinson and Company.



ELISA v.s. CBA

Comparison


ELISA	方法	BD CBA Flex Set	BD CBA Kit
1 個樣本 1 種	同時測定蛋白種類數目 勝	1 個樣本最多可 30 種	1 個樣本最多可 6 種
化學發光、吸光值；螢光	檢測原理	螢光	螢光
繁雜，需清洗 10 次以上。	實驗操作 勝	簡便，只需清洗一次。	
100 µl 以上	需要樣本量 勝	25-50 µl	
用 96-well plate	測定方式 勝	可用 96-well plate 或 Tube 來進行	
Plate Reader	測定儀器	2-Laser Flow Cytometry	1-Laser Flow Cytometry


Correlation

- *E. Morgan et al. / Clinical Immunology 110 (2004) 252–266*
- *BLOOD, 15 AUGUST 2005 VOLUME 106, NUMBER 4*
- *Electrophoresis 2009, 30, 4008–4019*
- *Methods, Volume 38, Issue 4, April 2006, Pages 317–323*
- *Cytometry Part B (Clinical Cytometry) 61B:35–39 (2004)*
- *Cytokine, Volume 54, Issue 2, May 2011, Pages 136–143*

..... etc

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.





BD CBA : designed for easy and efficient multiplexing

- Analyzed multiple analytes simultaneously.
- Reduced sample volume requirement.
- Reduced hands-on time by parallel analysis of samples.
- Wide dynamic range of fluorescence detection. (requires fewer sample dilutions)
- Offer automated sample acquisition and increased throughput.

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Beads Provide a Flexible Platform




**Different fluorescence
intensities**

Different colors with different intensities

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



How A Bead-based Immunoassay Works



© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Analytes

- Inflammatory mediators
- Chemokines
- Immunoglobulin isotypes
- Intracellular signaling molecules
- Apoptotic mediators
- Adhesion molecules
- Antibodies
-etc.

Ref.

Am J Clin Pathol 2002;118:346-353

Curr Protoc Cytom. 2006 Feb;Chapter 12:Unit 13.

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company. *Clinical Immunology* Volume 110, Issue 3, March 2004, Pages 272-276




一根釣竿就可以釣起所有想要的魚

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company. 

BD™ Cytometric Bead Array (CBA)


- Flow cytometry application
- Bead-Based Immunoassays
- Multiplexed quantitation and detection
- Up to 30 proteins can be analyzed using just 25 to 50 μL of sample.

BD™ CBA Flex Set を用いたタンパクレベル解析
 ヒト末梢血単核細胞 (PBMC) を異なる2つの条件下で刺激し、上清を回収後BD™ CBA Flex Set assay (30 plex) で測定した。



© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.


Assay Principle



© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Assay Procedure



© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.





Product Format

BD CBA Flex Sets

提供開放式的彈性實驗設計，可自由搭配想要偵測的試劑。

BD CBA Kits

將常常一起測試的 CBA 試劑合併成套裝組合，例如，Th1/Th2 cytokine kit、Inflammation Kit，可節省時間。


© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



CBA Kit

將常常同時偵測的種類，合併成套裝組合


以 Human BD CBA Th1/Th2 Cytokine Kit 實驗為例說明實驗流程 (用 BD Calibur 分析)。



© 2016 BD. BD, the


Capture beads 抓到想要測定的 cytokine 後，再跟帶有 PE 螢光抗體的 Detection reagent 以 Sandwich 方式作用。當測定的 cytokine 濃度越高，FL2 的螢光強度越強，Beads 群落會越往右偏移。





所得結果再以 FCAP Array 軟體分析之後，標準品曲線如上。

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.





Product Format

BD CBA Flex Sets	提供開放式的彈性實驗設計，可自由搭配想要偵測的試劑。
BD CBA Kits	將常會一起測試的 CBA 試劑合併成套裝組合，例如，Th1/Th2 cytokine kit、Inflammation Kit，可節省時間。

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



CBA FLEX SET



© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




CBA FLEX SET Product Line

- **Soluble Protein Flex Set**
 - Human Soluble Protein Flex Set
 - Mouse Soluble Protein Flex Set
 - Rat Soluble Protein Flex Set
- **Human TGF- β 1 Single Plex Flex Set**
- **Human Immunoglobulin Flex Set** (*For Ab Isotyping*)
- **Enhanced Sensitivity Flex Set** (*0.274–200 pg/mL*)
 - Human Enhanced Sensitivity Flex Set
 - Mouse Enhanced Sensitivity Flex Set
- **Cell Signaling Flex Set**
- **Functional Beads**

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Enhanced Sensitivity Flex Set System




Features

- Quantitation range of **0.274–200 pg/mL**
- Designed for multiplex analysis, providing more data using a single sample
- Open and configurable to maximize flexibility
- Analyzes soluble proteins with just 25 to 50 μL of sample
- Works with flow cytometers equipped with 488-nm and 633-nm lasers

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Time course of human PBMCs stimulated with PMA and Calcium Ionophore




© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Cell Signaling Flex Set System

- 定性及定量
 - 可同時偵測Cell lysate 中多種磷酸化蛋白
 - 標準品
- 單一、少量的樣本可進行多種分析
- 樣本製備與Western Blot 類似



Analysis of phospho-ERK1/2 protein levels in HEK 293 and NIH 3T3 cells in response to EGF, TNF, and PDGF stimulation.

A shows the results of Western blot analysis.

B shows the results of BD CBA Flex Set analysis.

Data courtesy of Dr. Tony Pawson and Dr. Jay Park, Mount Sinai Hospital, Toronto, Canada.

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Cell Signaling Flex Set System




Figure 2. Jurkat cells activated with anti-CD3 and anti-CD28 and crosslinked by anti-mouse Ig for 2 minutes. The reaction was stopped by the addition of SDS (1% final) whereupon the samples were placed in a boiling water bath. A phosphorylated recombinant standard was used to generate a standard curve. The upper panel shows the standard curve versus a titration of the activated Jurkat lysate. An important characteristic of the assay is that these two curves are parallel to each other. This ensures linearity of the readings so that twice as much lysate will have twice as many units/ml. The same concentrations of lysate were also run on an SDS-PAGE gel followed by immunoblotting with an anti-phospho-ZAP-70 antibody. Both a 30 second and a 10 minute exposure are shown. The BD CBA Flex Set assay is at least as sensitive as a Western blot. In this particular example, phosphorylation of ZAP-70 can be detected in a BD CBA Flex Set assay using lysate from less than 1,000 cells.

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




Human soluble protein 30 plex

Bead Position	4	5	6	7	8	9
A	IL-2	IL-4	IL-5	IL-6	IL-7	IL-8
B	IL-1 β	IP-10	IL-9	IL-10	VEGF	MIP-1 α
C	Angiogenin basic FGF		FasL	Eotaxin	G-CSF	GM-CSF
D	RANTES	IL-3	OSM	LT- α	MCP-1	TNF
E	MIP-1 β	IL-12p70	IL-13	IFN- γ	MIG	IgE

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.




Standard Curve



© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



BD™ CBA Functional Bead




© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company. DTT, Dithiothreitol; Sulfo-SMCC, Sulfosuccinimidyl 4-N-maleimidomethyl cyclohexane-1-carboxylate; NEM, N-Ethylmaleimide.



BD™ CBA Functional Bead

実験手順 (約4時間)



- Step 1: Bead Preparation
- Step 2: Protein Modification
- Step 3: Buffer Exchange to Remove Unreacted Components
- Step 4: Protein Conjugation

© 2016

company.



Advantages

- Simple conjugation procedure completed in less than 4 hours.
- Requires less than 100 µg of protein at a concentration of 1 mg/mL.
- Ability to conjugate any protein molecule containing a free amino group to the beads.
- Compatible with a wide selection of flow cytometers for ease of analysis.
- Specific reagents available for confirming the success of conjugation reactions.
- Supporting reagents and procedures available for performing instrument setups and assays

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Thank you!



90 © 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.