

*審查資料區間:2022/01-2022/12

允收標準參考出處:

A:原廠

B:醫檢學會定量檢驗品管作業指引裡的CLIA 規範每項分析量測的精密性的表現要求

C:醫檢學會定量檢驗品管作業指引裡的檢驗分析物的生理變異之個體內生理變異(CVw)

品管代號	組別	機台	項目代號	項目說明	Level	Mean	USD	aSD	Avg	CV%(aSD/Avg)	%	CV%(aSD/Mean)	原廠CV	允收標準參考出處	允收標準	允收標準	合格
5A	桃園組	H8AA	72-303	ALB	1	2.6	0.04	0.07	2.58	2.53	1.15%	2.513278347	2.00%	B	8.00%	B/2	✓
5A	桃園組	H8AA	72-303	ALB	3	4.4	0.06	0.05	4.33	1.06	0.85%	1.045789789	2.00%	B	8.00%	B/2	✓
5A	桃園組	H8AA	72-304	LDL	1	65	1.5	1.09	64.92	1.68	1.68%	1.673868192	5.00%	B	20.00%	B/2	✓
5A	桃園組	H8AA	72-304	LDL	3	163	3	3.23	161.64	2.00	2.00%	1.980390457	5.00%	B	20.00%	B/2	✓
5A	桃園組	H8AA	72-305	Lp:TCHO	1	163	2	1.50	168.38	0.89	0.89%	0.918786312	3.00%	B	10.00%	B/2	✓
5A	桃園組	H8AA	72-305	Lp:TCHO	2	316	4.5	3.23	328.48	0.98	0.98%	1.020938934	3.00%	B	10.00%	B/2	✓
5A	桃園組	H8AA	72-307	BUN	1	14.4	0.5	0.29	14.52	1.98	1.98%	1.997736587	5.00%	B	9.00%	B/2	✓
5A	桃園組	H8AA	72-307	BUN	3	69	1	1.16	67.23	1.73	1.73%	1.68632526	5.00%	B	9.00%	B/2	✓
5A	桃園組	H8AA	72-310	TG	1	101	2	1.47	102.98	1.43	1.43%	1.453302746	3.00%	B	15.00%	B/2	✓
5A	桃園組	H8AA	72-310	TG	3	241	4	6.54	232.84	2.81	2.81%	2.713044884	3.00%	B	15.00%	B/2	✓
5A	桃園組	H8AA	72-314	GLU	1	60	1.5	1.47	56.88	2.58	2.58%	2.449132382	5.00%	B	8.00%	B/2	✓
5A	桃園組	H8AA	72-314	GLU	3	357	5	6.60	344.91	1.91	1.91%	1.849430857	5.00%	B	8.00%	B/2	✓
5A	桃園組	H8AA	72-325	TP	1	4.06	0.1	0.12	3.77	3.08	3.08%	2.854266597	5.00%	B	8.00%	B/2	✓
5A	桃園組	H8AA	72-325	TP	3	6.7	0.1	0.15	6.60	2.22	2.22%	2.187721458	5.00%	B	8.00%	B/2	✓
5A	桃園組	H8AA	72-327	CA	1	5.73	0.15	0.11	5.74	1.90	1.82%	1.898160314	2.00%	B	+/- 1.0 mg/dL	B/2	✓
5A	桃園組	H8AA	72-327	CA	3	13.4	0.2	0.29	12.96	2.21	1.47%	2.138745568	2.00%	B	+/- 1.0 mg/dL	B/2	✓
5A	桃園組	H8AA	72-329	P	1	1.81	0.08	0.06	1.82	3.07	3.03%	3.079092627		B	10.00%	B/2	✓
5A	桃園組	H8AA	72-329	P	3	7.19	0.1	0.13	7.03	1.82	1.50%	1.776668992	1.75%	B	10.00%	B/2	✓
5A	桃園組	H8AA	72-331	UA	1	3.42	0.1	0.13	3.16	4.22	4.22%	3.906252994	7.00%	B	10.00%	B/2	✓
5A	桃園組	H8AA	72-331	UA	3	9.7	0.1	0.12	9.61	1.29	1.29%	1.28237794	7.00%	B	10.00%	B/2	✓
5A	桃園組	H8AA	72-333	CRE	1	0.77	0.04	0.03	0.75	3.98	3.38%	3.87431263		B	10.00%	B/2	✓
5A	桃園組	H8AA	72-333	CRE	3	6.5	0.1	0.12	6.51	1.88	1.70%	1.880075757	3.19%	B	10.00%	B/2	✓
5A	桃園組	H8AA	72-355	NA	1	117.7	1	1.71	115.37	1.48	0.83%	1.455014963	1.00%	A	1.00%	A	✓
5A	桃園組	H8AA	72-355	NA	3	158	1	0.96	158.91	0.60	0.59%	0.606825022	1.00%	A	1.00%	A	✓
5A	桃園組	H8AA	72-356	K	1	2.4	0.07	0.07	2.55	2.72	0.98%	2.889282354	2.00%	A	2.00%	A	✓
5A	桃園組	H8AA	72-356	K	3	7.54	0.1	0.09	7.54	1.23	1.45%	1.232545427	2.00%	A	2.00%	A	✓
5A	桃園組	H8AA	72-357	CL	1	76.2	1	1.28	74.00	1.73	1.73%	1.682189397	4.00%	A	4.00%	A/2	✓
5A	桃園組	H8AA	72-357	CL	3	121	1	1.05	120.90	0.87	0.87%	0.866911453	4.00%	A	4.00%	A/2	✓
5A	桃園組	H8AA	72-360	AST	1	38	1.5	1.62	34.16	4.75	4.75%	4.268525421	10.00%	C	12.30%	C/2	✓
5A	桃園組	H8AA	72-360	AST	3	212	3	3.39	207.82	1.63	1.63%	1.598713023	10.00%	C	12.30%	C/2	✓
5A	桃園組	H8AA	72-361	ALT	1	26	1.4	1.43	28.74	4.99	4.99%	5.511204563	10.00%	B	15.00%	B/2	✓
5A	桃園組	H8AA	72-361	ALT	3	184	2.5	2.35	182.41	1.29	1.29%	1.275267421	10.00%	B	15.00%	B/2	✓
5A	桃園組	H8AA	72-363	ALP	1	25	2	1.90	23.83	7.99	6.83%	7.615508394		B	20.00%	B/2	✓
5A	桃園組	H8AA	72-363	ALP	3	274	7	11.72	260.00	4.51	2.20%	4.276494227	3.85%	B	20.00%	B/2	✓
5A	桃園組	H8AA	72-365	T-BIL	1	0.54	0.05	0.02	0.60	3.84	3.84%	4.274479875	5.00%	B	20.00%	B/2	✓
5A	桃園組	H8AA	72-365	T-BIL	3	6.8	0.1	0.12	6.83	1.70	1.70%	1.709500433	5.00%	B	20.00%	B/2	✓
5A	桃園組	H8AA	72-366	D-BIL	1	0.29	0.04	0.01	0.28	4.57	4.57%	4.416408644	5.00%	C	36.80%	C/2	✓
5A	桃園組	H8AA	72-366	D-BIL	3	1.58	0.05	0.04	1.65	2.70	2.70%	2.827356716	5.00%	C	36.80%	C/2	✓
5A	桃園組	H8AA	72-368	rGT	1	23.5	1	0.60	23.11	2.61	2.61%	2.570822249	5.00%	B	15.00%	B/2	✓
5A	桃園組	H8AA	72-368	rGT	3	114	2	2.40	111.85	2.15	2.15%	2.10529333	5.00%	B	15.00%	B/2	✓
5A	桃園組	H8AA	72-375	Homocyst	1	13.1	0.4	0.44	12.80	3.43	2.90%	3.353677584	2.91%	C	8.30%	C/2	✓
5A	桃園組	H8AA	72-375	Homocyst	2	33.1	1.6	1.54	32.96	4.67	3.03%	4.650282242		C	8.30%	C/2	✓
5A	桃園組	H8AA	72-434	Lp(a)	1	9.5	0.3	0.28	10.18	2.72	2.72%	2.912359314	4.76%	C	20.80%	C/2	✓
5A	桃園組	H8AA	72-434	Lp(a)	2	18.2	0.6	0.39	20.04	1.93	1.93%	2.130615296	4.76%	C	20.80%	C/2	✓
5A	桃園組	H8AA	72-547	CRP	1	0.83	0.1	0.05	1.18	4.34	4.34%	6.154849637	10.00%	C	42.20%	C/2	✓
5A	桃園組	H8AA	72-547	CRP	2	5.11	0.15	0.13	4.76	2.69	2.69%	2.502279398	10.00%	C	42.20%	C/2	✓
5A	桃園組	H8AA	72A428	HDL-C	1	27	1	0.42	24.20	1.72	1.72%	1.538819184	5.00%	B	20.00%	B/2	✓
5A	桃園組	H8AA	72A428	HDL-C	2	57	2	1.05	52.44	2.01	2.01%	1.848157942	5.00%	B	20.00%	B/2	✓
5B	桃園組	H8BA	72-303	ALB	1	2.6	0.04	0.08	2.56	2.96	1.07%	2.919105502	2.00%	B	8.00%	B/2	✓
5B	桃園組	H8BA	72-303	ALB	3	4.4	0.06	0.04	4.34	1.01	0.97%	0.991370063	2.00%	B	8.00%	B/2	✓
5B	桃園組	H8BA	72-304	LDL	1	65	1.5	1.20	65.98	1.81	1.81%	1.841082942	5.00%	B	20.00%	B/2	✓
5B	桃園組	H8BA	72-304	LDL	3	163	3	3.24	163.53	1.98	1.98%	1.98952455	5.00%	B	20.00%	B/2	✓
5B	桃園組	H8BA	72-305	Lp:TCHO	1	163	2	1.63	167.83	0.97	0.97%	1.003061103	3.00%	B	10.00%	B/2	✓

*審查資料區間:2022/01-2022/12

允收標準參考出處:

A:原廠

B:醫檢學會定量檢驗品管作業指引裡的CLIA 規範每項分析量測的精密性的表現要求

C:醫檢學會定量檢驗品管作業指引裡的檢驗分析物的生理變異之個體內生理變異(CVw)

品管代號	組別	機台	項目代號	項目說明	Level	Mean	USD	aSD	Avg	CV%(aSD/Avg)	%	CV%(aSD/Mean)	原廠CV	允收標準參考出處	允收標準	允收標準	合格
5B	桃園組	H8BA	72-305	Lp:TCHO	2	316	4.5	3.45	327.55	1.05	1.05%	1.090877579	3.00%	B	10.00%	B/2	✓
5B	桃園組	H8BA	72-307	BUN	1	14.4	0.5	0.30	14.56	2.05	2.05%	2.077102911	5.00%	B	9.00%	B/2	✓
5B	桃園組	H8BA	72-307	BUN	3	69	1	1.17	67.17	1.74	1.74%	1.696546768	5.00%	B	9.00%	B/2	✓
5B	桃園組	H8BA	72-310	TG	1	101	2	1.54	104.37	1.47	1.47%	1.522690801	3.00%	B	15.00%	B/2	✓
5B	桃園組	H8BA	72-310	TG	3	241	4	6.37	232.82	2.73	2.73%	2.641462999	3.00%	B	15.00%	B/2	✓
5B	桃園組	H8BA	72-314	GLU	1	60	1.5	1.25	57.63	2.17	2.17%	2.085515377	5.00%	B	8.00%	B/2	✓
5B	桃園組	H8BA	72-314	GLU	3	357	5	6.12	345.94	1.77	1.77%	1.71463135	5.00%	B	8.00%	B/2	✓
5B	桃園組	H8BA	72-325	TP	1	4.06	0.1	0.12	3.77	3.27	3.27%	3.035474438	5.00%	B	8.00%	B/2	✓
5B	桃園組	H8BA	72-325	TP	3	6.7	0.1	0.15	6.62	2.22	2.22%	2.189571559	5.00%	B	8.00%	B/2	✓
5B	桃園組	H8BA	72-327	CA	1	5.73	0.15	0.10	5.78	1.74	1.79%	1.761000775	2.00%	B	+/- 1.0 mg/dL	B/2	✓
5B	桃園組	H8BA	72-327	CA	3	13.4	0.2	0.27	12.96	2.07	1.44%	2.002577995	2.00%	B	+/- 1.0 mg/dL	B/2	✓
5B	桃園組	H8BA	72-329	P	1	1.81	0.08	0.06	1.83	3.19	3.03%	3.229868831		B	10.00%	B/2	✓
5B	桃園組	H8BA	72-329	P	3	7.19	0.1	0.12	7.09	1.63	1.51%	1.606554441	1.75%	B	10.00%	B/2	✓
5B	桃園組	H8BA	72-331	UA	1	3.42	0.1	0.11	3.16	3.60	3.60%	3.329667331	7.00%	B	10.00%	B/2	✓
5B	桃園組	H8BA	72-331	UA	3	9.7	0.1	0.11	9.72	1.17	1.17%	1.169333407	7.00%	B	10.00%	B/2	✓
5B	桃園組	H8BA	72-333	CRE	1	0.77	0.04	0.03	0.76	3.83	3.40%	3.801222262		B	10.00%	B/2	✓
5B	桃園組	H8BA	72-333	CRE	3	6.5	0.1	0.11	6.51	1.72	1.59%	1.72784664	3.19%	B	10.00%	B/2	✓
5B	桃園組	H8BA	72-355	NA	1	117.7	1	1.86	115.62	1.61	0.81%	1.578229788	1.00%	A	1.00%	A	✓
5B	桃園組	H8BA	72-355	NA	3	158	1	0.93	158.34	0.59	0.55%	0.589607367	1.00%	A	1.00%	A	✓
5B	桃園組	H8BA	72-356	K	1	2.4	0.07	0.06	2.57	2.30	1.56%	2.456112421	2.00%	A	2.00%	A	✓
5B	桃園組	H8BA	72-356	K	3	7.54	0.1	0.09	7.52	1.16	0.87%	1.157657802	2.00%	A	2.00%	A	✓
5B	桃園組	H8BA	72-357	CL	1	76.2	1	1.33	73.81	1.81	1.81%	1.751219452	4.00%	A	4.00%	A/2	✓
5B	桃園組	H8BA	72-357	CL	3	121	1	0.97	120.47	0.81	0.81%	0.804028857	4.00%	A	4.00%	A/2	✓
5B	桃園組	H8BA	72-360	AST	1	38	1.5	1.61	34.80	4.62	4.62%	4.226787156	10.00%	C	12.30%	C/2	✓
5B	桃園組	H8BA	72-360	AST	3	212	3	3.83	210.44	1.82	1.82%	1.808061038	10.00%	C	12.30%	C/2	✓
5B	桃園組	H8BA	72-361	ALT	1	26	1.4	1.36	29.16	4.66	4.66%	5.229060502	10.00%	B	15.00%	B/2	✓
5B	桃園組	H8BA	72-361	ALT	3	184	2.5	2.63	183.52	1.43	1.43%	1.430065352	10.00%	B	15.00%	B/2	✓
5B	桃園組	H8BA	72-363	ALP	1	25	2	1.70	24.75	6.88	6.23%	6.813729741		B	20.00%	B/2	✓
5B	桃園組	H8BA	72-363	ALP	3	274	7	12.12	260.37	4.65	1.94%	4.423084179	3.85%	B	20.00%	B/2	✓
5B	桃園組	H8BA	72-365	T-BIL	1	0.54	0.05	0.03	0.60	4.17	4.17%	4.641949634	5.00%	B	20.00%	B/2	✓
5B	桃園組	H8BA	72-365	T-BIL	3	6.8	0.1	0.13	6.78	1.90	1.90%	1.896115919	5.00%	B	20.00%	B/2	✓
5B	桃園組	H8BA	72-366	D-BIL	1	0.29	0.04	0.01	0.28	4.75	4.75%	4.601087363	5.00%	C	36.80%	C/2	✓
5B	桃園組	H8BA	72-366	D-BIL	3	1.58	0.05	0.05	1.66	2.93	2.93%	3.084699204	5.00%	C	36.80%	C/2	✓
5B	桃園組	H8BA	72-368	rGT	1	23.5	1	0.59	21.92	2.69	2.69%	2.507297419	5.00%	B	15.00%	B/2	✓
5B	桃園組	H8BA	72-368	rGT	3	114	3	3.72	108.43	3.43	3.43%	3.258856409	5.00%	B	15.00%	B/2	✓
5B	桃園組	H8BA	72-547	CRP	1	0.83	0.1	0.06	1.17	5.30	5.30%	7.476158114	10.00%	C	42.20%	C/2	✓
5B	桃園組	H8BA	72-547	CRP	2	5.11	0.15	0.13	4.82	2.62	2.62%	2.467423578	10.00%	C	42.20%	C/2	✓
5B	桃園組	H8BA	72A428	Lp:HDL	1	27	1	0.45	24.25	1.85	1.85%	1.660483858	4.76%	B	20.00%	B/2	✓
5B	桃園組	H8BA	72A428	Lp:HDL	2	57	2	1.18	52.75	2.23	2.23%	2.065424856	4.76%	B	20.00%	B/2	✓
5I	桃園組	Roche-H(auto)	72-135	CA199	1	17.1	1.26	0.92	18.56	4.98	4.98%	5.40078972	10.00%	C	15.95%	C/2	✓
5I	桃園組	Roche-H(auto)	72-135	CA199	2	94.3	4.67	3.27	91.84	3.56	3.56%	3.467479232	10.00%	C	15.95%	C/2	✓
5I	桃園組	Roche-H(auto)	72-138	CEA	1	4.63	0.15	0.36	4.38	8.12	2.09%	7.684610998	5.00%	B	15.00%	B/2	✓
5I	桃園組	Roche-H(auto)	72-138	CEA	2	48.5	2.59	1.64	49.53	3.31	1.83%	3.375797249	5.00%	B	15.00%	B/2	✓
5I	桃園組	Roche-H(auto)	72-142	HBSAg	1	0.38	0.08	0.04	0.38	9.73	9.73%	9.719980009	10.00%	B	REACTIVE (POS) OR NONREACTIVE(NEG)	B	✓
5I	桃園組	Roche-H(auto)	72-142	HBSAg	2	3.82	0.202	0.32	3.81	8.38	8.38%	8.358616986	10.00%	B	REACTIVE (POS) OR NONREACTIVE(NEG)	B	✓
5I	桃園組	Roche-H(auto)	72-143	Anti-HBS	1	2	0.484	0.00	2.00	0.00	0.00%	0	8.00%	B	REACTIVE (POS) OR NONREACTIVE(NEG)	B	✓
5I	桃園組	Roche-H(auto)	72-143	Anti-HBS	2	102	3.4	3.46	92.40	3.74	3.74%	3.390784214	8.00%	B	REACTIVE (POS) OR NONREACTIVE(NEG)	B	✓
5I	桃園組	Roche-H(auto)	72-192	CYFRA21-1	1	2.7	0.15	0.15	2.74	5.35	5.35%	5.4368812	7.00%	C	22.20%	C/2	✓
5I	桃園組	Roche-H(auto)	72-192	CYFRA21-1	2	26.8	1.43	1.10	26.89	4.00	4.00%	4.116092311	7.00%	C	22.20%	C/2	✓
5I	桃園組	Roche-H(auto)	72-955	Anti-HCV	1	0.047	0.07	0	0.05	9.00	4.47%	9.956527242	8.00%	B	REACTIVE (POS) OR NONREACTIVE(NEG)	B	✓
5I	桃園組	Roche-H(auto)	72-955	Anti-HCV	2	3.37	0.26	0.28	3.91	7.20	2.88%	8.354956299	8.00%	B	REACTIVE (POS) OR NONREACTIVE(NEG)	B	✓
5J	桃園組	Alinity-H(auto)	72-137	CA153	2	38.9	2.04	2.46	39.39	6.23	4.63%	6.312412814	8.00%	C	6.10%	C	✓
5J	桃園組	Alinity-H(auto)	72-137	CA153	3	68.5	3.8	6.43	69.59	9.24	3.94%	9.385738355	8.00%	C	6.10%	C	✓

*審查資料區間:2022/01-2022/12

允收標準參考出處:

A:原廠

B:醫檢學會定量檢驗品管作業指引裡的CLIA 規範每項分析量測的精密性的表現要求

C:醫檢學會定量檢驗品管作業指引裡的檢驗分析物的生理變異之個體內生理變異(CVw)

品管代號	組別	機台	項目代號	項目說明	Level	Mean	USD	aSD	Avg	CV%(aSD/Avg)	%	CV%(aSD/Mean)	原廠CV	允收標準參考出處	允收標準	允收標準	合格
5J	桃園組	Alinity-H(auto)	72-139	PSA	2	3.9	0.2	0.16	3.82	4.25	4.25%	4.159895036	8.00%	B	20.00%	B/2	✓
5J	桃園組	Alinity-H(auto)	72-139	PSA	3	23	1.22	1.33	22.35	5.95	5.95%	5.7857121	8.00%	B	20.00%	B/2	✓
5J	桃園組	Alinity-H(auto)	72-140	AFP	1	5.4	0.19	0.26	5.68	4.60	4.60%	4.843642097	7.50%	B	20.00%	B/2	✓
5J	桃園組	Alinity-H(auto)	72-140	AFP	3	200	9.82	9.76	196.22	4.98	4.98%	4.881699349	7.50%	B	20.00%	B/2	✓
5J	桃園組	Alinity-H(auto)	72-141	SCC	1	2	0.13	0.08	2.02	4.02	4.02%	4.071243513	10.00%	C	39.40%	C/2	✓
5J	桃園組	Alinity-H(auto)	72-141	SCC	2	50	3.3	1.82	51.49	3.53	3.53%	3.633176057	10.00%	C	39.40%	C/2	✓
5J	桃園組	Alinity-H(auto)	72-193	CA125	2	35	1.72	1.58	35.21	4.49	4.49%	4.518918291	10.00%	B	20.00%	B/2	✓
5J	桃園組	Alinity-H(auto)	72-193	CA125	3	78.6	3.46	3.61	74.70	4.83	4.83%	4.593885536	10.00%	B	20.00%	B/2	✓
5J	桃園組	Alinity-H(auto)	72-389	Total T4	2	6.6	0.34	0.30	6.79	4.40	4.40%	4.522298978	10.00%	A	10.00%	A	✓
5J	桃園組	Alinity-H(auto)	72-389	Total T4	3	13	0.68	1.14	13.48	8.45	8.45%	8.757348097	10.00%	A	10.00%	A	✓
5J	桃園組	Alinity-H(auto)	72-390	Free T4	1	0.77	0.04	0.06	0.74	7.96	7.96%	7.691754983	10.00%	A	10.00%	A	✓
5J	桃園組	Alinity-H(auto)	72-390	Free T4	3	2.58	0.17	0.17	2.61	6.54	6.54%	6.611977013	10.00%	A	10.00%	A	✓
5J	桃園組	Alinity-H(auto)	72-391	T3	1	81	3.7	3.06	82.81	3.70	3.70%	3.777724632	10.00%	C	9.40%	C/2	✓
5J	桃園組	Alinity-H(auto)	72-391	T3	3	279	13.2	11.28	271.26	4.16	4.16%	4.044614016	10.00%	C	9.40%	C/2	✓
5J	桃園組	Alinity-H(auto)	72-393	TSH	1	0.058	0.004	0.01	0.05	9.93	9.93%	9.136241078	10.00%	C	19.30%	C	✓
5J	桃園組	Alinity-H(auto)	72-393	TSH	3	19.2	0.97	0.78	19.25	4.05	4.05%	4.060784018	10.00%	C	19.30%	C	✓
W1	桃園組	Premier A1c1(H)	72-317	HbA1C	1	6	0.07	0.08	6.00	1.30	1.30%	1.301352461	3.00%	B	10.00%	B/2	✓
W1	桃園組	Premier A1c1(H)	72-317	HbA1C	2	10.65	0.18	0.17	10.51	1.61	1.61%	1.58784779	3.00%	B	10.00%	B/2	✓
W3	桃園組	Premier A1c2(H)	72-317	HbA1C	1	6	0.07	0.08	6.00	1.37	1.37%	1.369331779	3.00%	B	10.00%	B/2	✓
W3	桃園組	Premier A1c2(H)	72-317	HbA1C	2	10.65	0.18	0.16	10.53	1.56	1.56%	1.539319006	3.00%	B	10.00%	B/2	✓
W6	桃園組	CS-1600(桃園)	72A037	PT	1	11.6	0.25	0.25	11.77	2.14	2.14%	2.171333058	5.00%	B	15.00%	B/2	✓
W6	桃園組	CS-1600(桃園)	72A037	PT	3	61.1	2.5	1.87	68.13	2.74	2.74%	3.057256218	5.00%	B	15.00%	B/2	✓
W6	桃園組	CS-1600(桃園)	72A038	APTT	1	27.4	1	0.55	27.49	1.99	1.99%	1.997057374	5.00%	B	15.00%	B/2	✓
W6	桃園組	CS-1600(桃園)	72A038	APTT	3	58.9	1.6	1.01	64.64	1.56	1.56%	1.707949248	5.00%	B	15.00%	B/2	✓
W8	桃園組	ESR(SRS20II)(桃園)	72-013	ESR	1	5.5	1.5	1.43	5.62	25.39	25.39%	25.94543881	30.00%	A	30.00%	A	✓
W8	桃園組	ESR(SRS20II)(桃園)	72-013	ESR	2	60	4	1.99	60.49	3.29	3.29%	3.312841732	30.00%	A	30.00%	A	✓
WI	桃園組	XN-9000-A(桃園)	72A001	WBC	1	3.16	0.1	0.09	3.17	2.85	2.85%	2.857316634	10.00%	C	11.45%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72A001	WBC	2	7.07	0.16	0.21	7.26	2.96	2.96%	3.040305358	10.00%	C	11.45%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72A001	WBC	3	16.34	0.3	0.36	16.80	2.17	2.17%	2.228457613	10.00%	C	11.45%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72A015	Neut	1	39.7	1.2	1.27	38.61	3.30	3.30%	3.208113571	15.00%	C	17.10%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72A015	Neut	2	43.6	1.1	1.16	43.62	2.67	2.67%	2.668552054	15.00%	C	17.10%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72A015	Neut	3	47.4	1	0.96	47.79	2.01	2.01%	2.022794707	15.00%	C	17.10%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72B001	RBC	1	2.36	0.03	0.03	2.35	1.45	1.45%	1.439320551	5.00%	B	4.00%	B	✓
WI	桃園組	XN-9000-A(桃園)	72B001	RBC	2	4.28	0.04	0.08	4.27	1.94	1.94%	1.937600163	5.00%	B	4.00%	B	✓
WI	桃園組	XN-9000-A(桃園)	72B001	RBC	3	5.05	0.06	0.08	4.99	1.55	1.55%	1.528944059	5.00%	B	4.00%	B	✓
WI	桃園組	XN-9000-A(桃園)	72C001	HGB	1	5.39	0.07	0.11	5.44	2.07	2.07%	2.093462917	5.00%	B	4.00%	B	✓
WI	桃園組	XN-9000-A(桃園)	72C001	HGB	2	11.3	0.14	0.10	11.37	0.88	0.88%	0.884650633	5.00%	B	4.00%	B	✓
WI	桃園組	XN-9000-A(桃園)	72C001	HGB	3	14.7	0.2	0.32	15.05	2.15	2.15%	2.20423342	5.00%	B	4.00%	B	✓
WI	桃園組	XN-9000-A(桃園)	72C015	LYM	1	33.8	1.5	1.68	34.23	4.92	4.92%	4.983567599	15.00%	C	10.20%	C	✓
WI	桃園組	XN-9000-A(桃園)	72C015	LYM	2	29.4	1.2	1.37	28.39	4.81	4.81%	4.648385557	15.00%	C	10.20%	C	✓
WI	桃園組	XN-9000-A(桃園)	72C015	LYM	3	24.3	1	0.96	23.30	4.14	4.14%	3.965373748	15.00%	C	10.20%	C	✓
WI	桃園組	XN-9000-A(桃園)	72D001	HCT	1	16.3	0.3	0.45	16.54	2.72	2.72%	2.76176001	5.00%	B	4.00%	B	✓
WI	桃園組	XN-9000-A(桃園)	72D001	HCT	2	33.1	0.5	0.76	33.59	2.26	2.26%	2.297845786	5.00%	B	4.00%	B	✓
WI	桃園組	XN-9000-A(桃園)	72D001	HCT	3	42.5	0.6	0.93	43.45	2.15	2.15%	2.197589871	5.00%	B	4.00%	B	✓
WI	桃園組	XN-9000-A(桃園)	72D015	MON	1	12.2	1.5	1.52	13.15	11.58	11.58%	12.47651722	40.00%	C	17.50%	C	✓
WI	桃園組	XN-9000-A(桃園)	72D015	MON	2	11.8	1.1	1.50	12.71	11.76	11.76%	12.66990761	40.00%	C	17.50%	C	✓
WI	桃園組	XN-9000-A(桃園)	72D015	MON	3	12.3	1	0.85	12.66	6.75	6.75%	6.948861532	40.00%	C	17.50%	C	✓
WI	桃園組	XN-9000-A(桃園)	72E001	MCV	1	69.2	1	1.30	70.58	1.84	1.84%	1.879508996	5.00%	A	5.00%	A/2	✓
WI	桃園組	XN-9000-A(桃園)	72E001	MCV	2	77.3	0.8	0.80	77.58	1.03	1.03%	1.036102805	5.00%	A	5.00%	A/2	✓
WI	桃園組	XN-9000-A(桃園)	72E001	MCV	3	84.2	1	1.60	85.88	1.87	1.87%	1.905776605	5.00%	A	5.00%	A/2	✓
WI	桃園組	XN-9000-A(桃園)	72E015	EOS	1	9.5	0.8	0.73	9.25	7.85	7.85%	7.639352584	15.00%	C	21.00%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72E015	EOS	2	10.5	0.8	0.76	10.48	7.21	7.21%	7.200182304	15.00%	C	21.00%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72E015	EOS	3	11.2	1	0.81	11.41	7.08	7.08%	7.215490104	15.00%	C	21.00%	C/2	✓

*審查資料區間:2022/01-2022/12

允收標準參考出處:

A:原廠

B:醫檢學會定量檢驗品管作業指引裡的CLIA 規範每項分析量測的精密性的表現要求

C:醫檢學會定量檢驗品管作業指引裡的檢驗分析物的生理變異之個體內生理變異(CVw)

品管代號	組別	機台	項目代號	項目說明	Level	Mean	USD	aSD	Avg	CV%(aSD/Avg)	%	CV%(aSD/Mean)	原廠CV	允收標準參考出處	允收標準	允收標準	合格
WI	桃園組	XN-9000-A(桃園)	72F001	MCH	1	22.9	0.4	0.41	23.28	1.77	1.77%	1.796539613	10.00%	A	10.00%	A/2	✓
WI	桃園組	XN-9000-A(桃園)	72F001	MCH	2	26.3	0.4	0.27	26.45	1.02	1.02%	1.027338306	10.00%	A	10.00%	A/2	✓
WI	桃園組	XN-9000-A(桃園)	72F001	MCH	3	29.1	0.4	0.40	29.55	1.36	1.36%	1.378243222	10.00%	A	10.00%	A/2	✓
WI	桃園組	XN-9000-A(桃園)	72F015	BAS	1	4.8	0.2	0.16	4.78	3.25	3.25%	3.239986131	15.00%	C	28.00%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72F015	BAS	2	4.8	0.2	0.12	4.80	2.52	2.52%	2.517751695	15.00%	C	28.00%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72F015	BAS	3	4.8	0.2	0.11	4.84	2.19	2.19%	2.212277245	15.00%	C	28.00%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72G001	MCHC	1	33.1	0.5	0.43	32.98	1.31	1.31%	1.305311093	10.00%	A	10.00%	A/2	✓
WI	桃園組	XN-9000-A(桃園)	72G001	MCHC	2	34	0.5	0.51	34.32	1.50	1.50%	1.511753293	10.00%	A	10.00%	A/2	✓
WI	桃園組	XN-9000-A(桃園)	72G001	MCHC	3	34.5	0.5	0.45	35.25	1.26	1.26%	1.29136178	10.00%	A	10.00%	A/2	✓
WI	桃園組	XN-9000-A(桃園)	72H001	RDW-CV	1	19.6	0.2	0.27	19.33	1.39	1.39%	1.372619034	10.00%	C	3.50%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72H001	RDW-CV	2	18.1	0.2	0.18	17.98	0.98	0.98%	0.971727231	10.00%	C	3.50%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72H001	RDW-CV	3	16.9	0.15	0.11	17.15	0.66	0.66%	0.670656062	10.00%	C	3.50%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72I001	PLT	1	87	7	5.06	89.80	5.64	5.64%	5.818973373	10.00%	B	25.00%	B/2	✓
WI	桃園組	XN-9000-A(桃園)	72I001	PLT	2	251	10	9.52	248.45	3.83	3.83%	3.793136111	10.00%	B	25.00%	B/2	✓
WI	桃園組	XN-9000-A(桃園)	72I001	PLT	3	575	12	12.86	563.70	2.28	2.28%	2.236433164	10.00%	B	25.00%	B/2	✓
WI	桃園組	XN-9000-A(桃園)	72J001	RDW-SD	1	47.8	0.5	0.68	48.10	1.42	1.42%	1.426809887	10.00%	C	3.50%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72J001	RDW-SD	2	50.4	0.5	0.57	49.55	1.15	1.15%	1.132138451	10.00%	C	3.50%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72J001	RDW-SD	3	51.3	0.6	0.68	51.29	1.32	1.32%	1.320895151	10.00%	C	3.50%	C/2	✓
WI	桃園組	XN-9000-A(桃園)	72L015	NRBC%	1	6	2	0.55	4.74	11.50	11.50%	9.093382602	25.00%	A	25.00%	A/2	✓
WI	桃園組	XN-9000-A(桃園)	72L015	NRBC%	2	6	2	0.40	6.07	6.54	6.54%	6.615731997	25.00%	A	25.00%	A/2	✓
WI	桃園組	XN-9000-A(桃園)	72L015	NRBC%	3	6	2	0.23	6.63	3.41	3.41%	3.770284637	25.00%	A	25.00%	A/2	✓
WI	桃園組	XN-9000-A(桃園)	72NRBC	NRBC	1	0.15	0.03	0.02	0.15	10.50	10.50%	10.49781318	25.00%	A	25.00%	A/2	✓
WI	桃園組	XN-9000-A(桃園)	72NRBC	NRBC	2	0.43	0.03	0.02	0.44	5.60	5.60%	5.732308311	25.00%	A	25.00%	A/2	✓
WI	桃園組	XN-9000-A(桃園)	72NRBC	NRBC	3	1.08	0.05	0.04	1.11	3.58	3.58%	3.691840162	25.00%	A	25.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72-PLF	PLT-F	1	84	7	2.29	84.93	2.70	2.70%	2.729512749	10.00%	B	25.00%	B/2	✓
WJ	桃園組	XN-9000-B(桃園)	72-PLF	PLT-F	2	247	10	13.54	260.41	5.20	5.20%	5.482266601	10.00%	B	25.00%	B/2	✓
WJ	桃園組	XN-9000-B(桃園)	72-PLF	PLT-F	3	550	12	12.67	553.22	2.29	2.29%	2.303228201	10.00%	B	25.00%	B/2	✓
WJ	桃園組	XN-9000-B(桃園)	72A001	WBC	1	3.16	0.1	0.10	3.07	3.30	3.30%	3.204541598	10.00%	C	11.45%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72A001	WBC	2	7.07	0.16	0.19	7.04	2.69	2.69%	2.678215325	10.00%	C	11.45%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72A001	WBC	3	16.34	0.3	0.26	16.41	1.56	1.56%	1.565556344	10.00%	C	11.45%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72A015	Neut	1	39.7	1.2	1.23	38.22	3.23	3.23%	3.105951605	15.00%	C	17.10%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72A015	Neut	2	43.6	1.1	1.21	43.55	2.77	2.77%	2.771500941	15.00%	C	17.10%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72A015	Neut	3	47.4	1	1.05	47.77	2.20	2.20%	2.218616342	15.00%	C	17.10%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72B001	RBC	1	2.36	0.03	0.03	2.38	1.38	1.38%	1.387426321	5.00%	B	4.00%	B	✓
WJ	桃園組	XN-9000-B(桃園)	72B001	RBC	2	4.28	0.04	0.09	4.30	2.01	2.01%	2.024408037	5.00%	B	4.00%	B	✓
WJ	桃園組	XN-9000-B(桃園)	72B001	RBC	3	5.05	0.06	0.09	5.06	1.84	1.84%	1.840941165	5.00%	B	4.00%	B	✓
WJ	桃園組	XN-9000-B(桃園)	72C001	HGB	1	5.39	0.07	0.10	5.45	1.82	1.82%	1.844017904	5.00%	B	4.00%	B	✓
WJ	桃園組	XN-9000-B(桃園)	72C001	HGB	2	11.3	0.14	0.12	11.38	1.07	1.07%	1.075507626	5.00%	B	4.00%	B	✓
WJ	桃園組	XN-9000-B(桃園)	72C001	HGB	3	14.7	0.2	0.32	15.03	2.14	2.14%	2.186866704	5.00%	B	4.00%	B	✓
C	桃園組	XN-9000-B(桃園)	72C015	LYM	1	33.8	1.5	1.86	34.75	5.34	5.34%	5.493279697	15.00%	C	10.20%	C	✓
WJ	桃園組	XN-9000-B(桃園)	72C015	LYM	2	29.4	1.2	1.37	28.46	4.81	4.81%	4.652614803	15.00%	C	10.20%	C	✓
WJ	桃園組	XN-9000-B(桃園)	72C015	LYM	3	24.3	1	0.80	23.35	3.41	3.41%	3.27577739	15.00%	C	10.20%	C	✓
WJ	桃園組	XN-9000-B(桃園)	72D001	HCT	1	16.3	0.3	0.39	16.70	2.33	2.33%	2.382062691	5.00%	B	4.00%	B	✓
WJ	桃園組	XN-9000-B(桃園)	72D001	HCT	2	33.1	0.5	0.79	33.88	2.34	2.34%	2.396435381	5.00%	B	4.00%	B	✓
WJ	桃園組	XN-9000-B(桃園)	72D001	HCT	3	42.5	0.6	0.92	43.58	2.12	2.12%	2.174036555	5.00%	B	4.00%	B	✓
WJ	桃園組	XN-9000-B(桃園)	72D015	MON	1	12.2	1.5	1.52	12.80	11.89	11.89%	12.46914343	40.00%	C	17.50%	C	✓
WJ	桃園組	XN-9000-B(桃園)	72D015	MON	2	11.8	1.1	1.42	12.84	11.06	11.06%	12.03699022	40.00%	C	17.50%	C	✓
WJ	桃園組	XN-9000-B(桃園)	72D015	MON	3	12.3	1	0.73	12.69	5.74	5.74%	5.929019793	40.00%	C	17.50%	C	✓
WJ	桃園組	XN-9000-B(桃園)	72E001	MCV	1	69.2	1	1.29	70.52	1.83	1.83%	1.860311913	5.00%	A	5.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72E001	MCV	2	77.3	0.8	0.81	77.37	1.05	1.05%	1.047266065	5.00%	A	5.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72E001	MCV	3	84.2	1	1.86	85.81	2.17	2.17%	2.213920238	5.00%	A	5.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72E015	EOS	1	9.5	0.8	0.74	9.47	7.82	7.82%	7.795925745	15.00%	C	21.00%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72E015	EOS	2	10.5	0.8	0.77	10.35	7.43	7.43%	7.32763337	15.00%	C	21.00%	C/2	✓

*審查資料區間:2022/01-2022/12

允收標準參考出處:

A:原廠

B:醫檢學會定量檢驗品管作業指引裡的CLIA 規範每項分析量測的精密性的表現要求

C:醫檢學會定量檢驗品管作業指引裡的檢驗分析物的生理變異之個體內生理變異(CVw)

品管代號	組別	機台	項目代號	項目說明	Level	Mean	USD	aSD	Avg	CV%(aSD/Avg)	%	CV%(aSD/Mean)	原廠CV	允收標準參考出處	允收標準	允收標準	合格
WJ	桃園組	XN-9000-B(桃園)	72E015	EOS	3	11.2	1	0.81	11.38	7.12	7.12%	7.237540739	15.00%	C	21.00%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72F001	MCH	1	22.9	0.4	0.45	23.10	1.96	1.96%	1.980217109	10.00%	A	10.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72F001	MCH	2	26.3	0.4	0.35	26.17	1.32	1.32%	1.316118171	10.00%	A	10.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72F001	MCH	3	29.1	0.4	0.78	29.69	2.62	2.62%	2.668982875	10.00%	A	10.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72F015	BAS	1	4.8	0.2	0.17	4.76	3.57	3.57%	3.534466773	15.00%	C	28.00%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72F015	BAS	2	4.8	0.2	0.13	4.80	2.75	2.75%	2.748927866	15.00%	C	28.00%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72F015	BAS	3	4.8	0.2	0.12	4.81	2.56	2.56%	2.567759791	15.00%	C	28.00%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72G001	MCHC	1	33.1	0.5	0.45	32.73	1.39	1.39%	1.372424321	10.00%	A	10.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72G001	MCHC	2	34	0.5	0.61	34.08	1.79	1.79%	1.798267193	10.00%	A	10.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72G001	MCHC	3	34.5	0.5	0.51	34.94	1.45	1.45%	1.468705546	10.00%	A	10.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72H001	RDW-CV	1	19.6	0.2	0.25	19.31	1.29	1.29%	1.26805342	10.00%	C	3.50%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72H001	RDW-CV	2	18.1	0.2	0.17	17.93	0.96	0.96%	0.946643408	10.00%	C	3.50%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72H001	RDW-CV	3	16.9	0.15	0.11	17.07	0.62	0.62%	0.628291423	10.00%	C	3.50%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72I001	PLT	1	87	7	4.97	89.10	5.58	5.58%	5.709921446	10.00%	B	25.00%	B/2	✓
WJ	桃園組	XN-9000-B(桃園)	72I001	PLT	2	251	10	8.61	245.84	3.50	3.50%	3.431518453	10.00%	B	25.00%	B/2	✓
WJ	桃園組	XN-9000-B(桃園)	72I001	PLT	3	575	12	12.86	559.75	2.30	2.30%	2.235851159	10.00%	B	25.00%	B/2	✓
WJ	桃園組	XN-9000-B(桃園)	72J001	RDW-SD	1	47.8	0.5	0.82	47.80	1.72	1.72%	1.715496795	10.00%	C	3.50%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72J001	RDW-SD	2	50.4	0.5	0.43	49.29	0.86	0.86%	0.843268539	10.00%	C	3.50%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72J001	RDW-SD	3	51.3	0.6	0.52	50.81	1.02	1.02%	1.011262175	10.00%	C	3.50%	C/2	✓
WJ	桃園組	XN-9000-B(桃園)	72L015	NRBC%	1	6	2	0.54	4.78	11.32	11.32%	9.022270598	25.00%	A	25.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72L015	NRBC%	2	6	2	0.38	6.11	6.28	6.28%	6.398524585	25.00%	A	25.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72L015	NRBC%	3	6	2	0.22	6.62	3.29	3.29%	3.625142333	25.00%	A	25.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72NRBC	NRBC	1	0.15	0.03	0.01	0.15	10.16	10.16%	9.926706862	25.00%	A	25.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72NRBC	NRBC	2	0.43	0.03	0.02	0.43	5.58	5.58%	5.579250596	25.00%	A	25.00%	A/2	✓
WJ	桃園組	XN-9000-B(桃園)	72NRBC	NRBC	3	1.08	0.05	0.04	1.08	3.59	3.59%	3.609443742	25.00%	A	25.00%	A/2	✓
WG	桃園組	OC sensor PLEDIA	72-732	OB-EIA	1	149	9	7.13724883	144.79	4.93	4.93%		10.00%	A	10.00%	A	✓
WG	桃園組	OC sensor PLEDIA	72-732	OB-EIA	2	76	5	4.654987775	72.57	6.41	6.41%		10.00%	A	10.00%	A	✓
WG	桃園組	OC sensor PLEDIA	72-732	OB-EIA	3	0.001	3	0.549602133	0.06	9.16	9.16%		10.00%	A	10.00%	A	✓
	桃園組	Triage	72-572	Trop I	1	0.45	0.11	0.06	0.36	15.71	16%		12.00%	B	30.00%	B	✓
	桃園組	Triage	72-572	Trop I	2	15.3	2.3	1.33	12.44	11.00	11%		12.00%	B	30.00%	B	✓
	桃園組	Triage	72-572	CK-MB	1	5.3	1.2	0.59	4.13	14.00	14%		12.00%	B	25.00%	B	✓
	桃園組	Triage	72-573	CK-MB	2	40.5	9	3.50	31.07	11.00	11%		12.00%	B	25.00%	B	✓
	桃園組	Triage	72-574	Myoglobin	1	97.4	14.6	4.81	79.33	6.00	6%		12.00%	C	17.60%	C	✓
	桃園組	Triage	72-575	Myoglobin	2	261	57	20.65	191.25	11.00	11%		12.00%	C	17.60%	C	✓
	桃園組	iSTAT	72-510	Blood gas pH	1	7.038	0.016	0.005	7.04	0.07	0.07%		0.50%	B	0.04	B/2	✓
	桃園組	iSTAT	72-510	Blood gas pH	2	7.463	0.016	0.004	7.46	0.06	0.06%		0.50%	B	0.04	B/2	✓
	桃園組	iSTAT	72-510	Blood gas pH	3	7.65	0.016	0.004	7.64	0.06	0.06%		0.50%	B	0.04	B/2	✓
	桃園組	iSTAT	72-510	Blood gas pO2	1	88	5	2.353	88.58	2.66	2.66%		6.00%	B	15.00%	B/2	✓
	桃園組	iSTAT	72-510	Blood gas pO2	2	121	6	4.209	117.08	3.60	3.60%		6.00%	B	15.00%	B/2	✓
	桃園組	iSTAT	72-510	Blood gas pO2	3	150	8	6.566	143.25	4.58	4.58%		6.00%	B	15.00%	B/2	✓
	桃園組	iSTAT	72-510	Blood gas pCO2	1	54.9	2.5	1.913	57.36	3.33	3.33%		5.00%	B	8.00%	B	✓
	桃園組	iSTAT	72-510	Blood gas pCO2	2	29	2.3	1.440	31.20	4.62	4.62%		5.00%	B	8.00%	B	✓
	桃園組	iSTAT	72-510	Blood gas pCO2	3	23	2.1	0.423	24.51	1.73	1.73%		5.00%	B	8.00%	B	✓